




# PureSet™

## Product overview

Discover advanced simplicity with the brand-new PureSet Trays. Small practices or large hospitals can enjoy the ease of use, and benefit from workflow improvements thanks to the robust stainless steel\* trays designed for automated reprocessing.

### Uncompromised hygiene

PureSet Trays are designed for automated cleaning with few manual steps. Tray and instruments stay together throughout the automated reprocessing workflow.



### Secure retainment of instruments

Instruments are contained and secured using metal springs.



### Organized to save you time

Straightforward layout with easily identifiable components designed to reflect the surgical workflow. Instruments are organized according to drilling protocol.

\*PureSet Trays are made of stainless steel except for the PureSet plates and grommets.

# NobelActive® PureSet 300402

# NobelParallel™ CC PureSet 300403

## Included in the kit


PUR0200	NobelActive®/NobelParallel™ CC PureSet Tray
34584	NobelActive Manual Torque Wrench Surgical <sup>1</sup>
32180	Surgical Driver
32112	Direction Indicator Ø 2/Ø 2.4–2.8 mm x 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
29164	Drill Extension Shaft
29532	Implant/Prosthetic Organizer
29543	Implant Sleeve Holder
37887	NobelActive Radiographic Template
300565	NobelActive PureSet Wall Chart

## Included in the kit

PUR0200	NobelActive®/NobelParallel™ CC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical <sup>1</sup>
32112	Direction Indicator Ø 2/Ø 2.4–2.8 mm x 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
29164	Drill Extension Shaft
29532	Implant/Prosthetic Organizer
29543	Implant Sleeve Holder
37792	NobelParallel Radiographic Template
300566	NobelParallel CC PureSet Wall Chart

### Implant Drivers

28 mm	36773	36718	36720 37859
37 mm	36774	36719	36721 37860



### Implant Drivers

28 mm	36718	36720 37859
37 mm	36719	36721 37860



All-in-one: Store NobelActive and NobelParallel CC surgical instruments in one single tray.

## Not included in the kit (available for separate purchase)

### Bone Mills with Guide

Ø 4.0	37861 <sup>2</sup>		
Ø 4.4	37863		
Ø 5.2	37864	37866	
Ø 6.2	37867		
Ø 6.7	37869		



### Bone Mill Guides

37862 <sup>2</sup>	37865	37868	37870



### Drills

Precision	36118
Ø 1.5	31278 (7–15 mm) <sup>2</sup>
Ø 2.0	32296 (7–10 mm) 32297 (7–15 mm) 32299 (10–18 mm)
Ø 2.4/2.8	32260 (7–10 mm) 32261 (7–15 mm) 32262 (10–18 mm)
Ø 2.8/3.2	37873 (7–10 mm) 34638 (7–15 mm) 34639 (10–18 mm)
Ø 3.2/3.6	32263 (7–10 mm) <sup>3</sup> 32264 (7–15 mm) 32265 (10–18 mm)
Ø 3.8/4.2	32275 (7–10 mm) <sup>3</sup> 32276 (7–15 mm) 32277 (10–18 mm)
Ø 4.2/4.6	37874 (7–10 mm) 34582 (7–15 mm) 34583 (10–18 mm)
Ø 4.2/5.0	37875 (7–10 mm) 37876 (7–15 mm)



### Screw Taps NobelActive

Ø 3.0	36816		
Ø 3.5	36236		
Ø 4.3	36237		
Ø 5.0	36238		
Ø 5.5	37871 (7–10 mm) 37872 (11.5–15 mm)		



### Screw Taps NobelParallel CC

Ø 3.75	37990 (7–13 mm) 37991 (7–18 mm)	
Ø 4.3	37992 (7–13 mm) 37993 (7–18 mm)	
Ø 5.0	37994 (7–13 mm) 37995 (7–18 mm)	
Ø 5.5	37996 (7–10 mm) 37997 (7–15 mm)	



### Cortical Drills NobelParallel CC

Ø 3.75	38000	
Ø 4.3	38001	
Ø 5.0	38002	
Ø 5.5	38003	




1 Article to be disassembled for cleaning  
2 For NobelActive only  
3 For NobelParallel CC only

# NobelReplace® CC PureSet 300546

## Included in the kit

PUR0300	NobelReplace® CC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical <sup>1</sup>
29149	Screwdriver Manual Unigrip 28 mm
29164	Drill Extension Shaft
29532	Implant/Prosthetic Organizer
29543	Implant Sleeve Holder
2042	Irrigation Needle
37320	NobelReplace Radiographic Template
300567	NobelReplace® CC PureSet Wall Chart



## Drills

	NP	RP	WP
Ø 3.5	32075 (8 mm) 29367 (10 mm) 36113 (11.5 mm) 29368 (13 mm) 29369 (16 mm)		
Ø 4.3		32076 (8 mm) 29370 (10 mm) 36114 (11.5 mm) 29371 (13 mm) 29372 (16 mm)	
Ø 5.0			32077 (8 mm) 29373 (10 mm) 36115 (11.5 mm) 29374 (13 mm) 29375 (16 mm)



## Dense Bone Drills

	NP	RP	WP
Ø 3.5	29377 (13 mm) 29378 (16 mm)		
Ø 4.3		29380 (13 mm) 29381 (16 mm)	
Ø 5.0			29383 (13 mm) 29384 (16 mm)



## Screw Taps

	NP	RP	WP
	36717	32090	32091



## Direction Indicators

	NP	RP	WP
	32255	32256	32257



## Implant Drivers

	NP	RP
28 mm	36718	36720
37 mm	36719	36721



## Not included in the kit (available for separate purchase)

### Bone Mills with Guide

	NP	RP
Ø 4.4	37863	
Ø 5.2	37864	37866
Ø 6.2		37867




### Bone Mill Guides

	NP	RP
	37865	37868




### Drills

36117	Drill with Tip Tapered 2 mm	
36118	Precision Drill	

Sterilization accessories <sup>4</sup>	
JN295	Dental Container 312 x 190 x 92 mm
JK091	Reusable Filter for Dental Container 2/pkg
JF166	Paper Single Use Filter for Dental Container 100/pkg
JG687	Indicator Label with Adhesive Back 100/pkg
JG752	Tamperproof Lock Orange Steam Indicator 100/pkg



Spare parts	
PUR0201	NobelActive/NobelParallel CC PureSet Plate
PUR0301	NobelReplace CC PureSet Plate



<sup>4</sup> Products available in selected markets only. For more information, please contact your sales representative.



### Order online

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

[nobelbiocare.com/store](http://nobelbiocare.com/store)



### Order by phone

Call our customer service team, or contact your sales representative.



On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit [nobelbiocare.com/warranty](http://nobelbiocare.com/warranty)



[nobelbiocare.com/pureset](http://nobelbiocare.com/pureset)

85003 GB 1804 Printed in the EU © Nobel Biocare Services AG, 2018. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](http://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. Disclaimer: Some products may not be regulatory cleared/ released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed dentist. See Instructions for Use for full prescribing information, including indications, contraindications, warnings and precautions.